

Anti-Hand2 Polyclonal Antibody For Wb, If, Ihc, Elisa - P61296

Item Number: CM0026822

Introduction

Anti-HAND2 polyclonal antibody validated for WB, IF/ICC, IHC-P, and ELISA. Reacts with human, mouse, rat HAND2. Suitable for cardiac morphogenesis and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	HAND2 Rabbit pAb
Remarks/Alias	Hed, dHand, DHAND2, Thing2, bHLHa26, HAND2
Species	Human
GeneID (Human)	9464
GeneID	9464
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61296
Protein Weight	24 kDa
Shipping	Ice bag