

Anti-Hao1/Gox Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ujm8

Item Number: CM0006662



Introduction

This rabbit monoclonal antibody specifically detects human HA01/GOX (2-hydroxyacid oxidase 1), a key peroxisomal enzyme in glyoxylate metabolism. Validated in WB, IHC-P, and ELISA, it cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	HA01/GOX Rabbit mAb
Remarks/Alias	GO; GOX; GOX1; HAOX1; HA01/GOX
Species	Human
GeneID (Human)	54363
GeneID	54363
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HA01/GOX (NP_060015.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UJM8
Protein Weight	41kDa
Shipping	Ice bag