

Anti-Hao1/Gox Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Ujm8

Item Number: CM0026760



Introduction

Rabbit polyclonal antibody targeting human HA01/GOX (glycolate oxidase) for applications including Western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects human, mouse, and rat HA01. Suitable for metabolic and kidney stone research.

[Learn More](#)

Parameter	Value
Product Name	HA01/GOX Rabbit pAb
Remarks/Alias	GO; GOX; GOX1; HAOX1; HA01/GOX
Species	Human
GeneID (Human)	54363
GeneID	54363
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJM8
Protein Weight	41 kDa
Shipping	Ice bag