

Anti-Hap1 Rabbit Polyclonal Antibody For Wb, Elisa - P54257

Item Number: CM0020691

Introduction

Rabbit polyclonal antibody against human Huntingtin-associated protein 1 (HAP1), validated for WB and ELISA, cross-reactive with mouse and rat. Ideal for studying HAP1's roles in neuronal trafficking, autophagy, and Huntington's disease.

[Learn More](#)

Parameter	Description
Product Name	HAP1 Rabbit pAb
Remarks/Alias	HLP; HAP2; HIP5; hHLP1; HAP1
Species	Human
GeneID (Human)	9001
GeneID	9001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human HAP1 (NP_817084.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P54257
Protein Weight	76kDa
Shipping	Ice bag