

Anti-Hapln1 Monoclonal Antibody For Wb And Elisa - P10915

Item Number: CM0004258

Introduction

Recombinant rabbit monoclonal antibody targeting human HAPLN1, a 40 kDa cartilage link protein. Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	HAPLN1 Rabbit mAb
Remarks/Alias	CRT1; CRTL1; HAPLN1
Species	Human
Gene ID (Human)	1404
Gene ID	1404
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-354 of human HAPLN1 (P10915)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10915
Protein Weight	40kDa
Shipping	Ice bag