


Anti-Haptoglobin (Hp) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P00738

Item Number: CM0004989



Introduction

High-quality rabbit monoclonal antibody against human Haptoglobin (HP), validated for WB, IF-P, IHC-P, and ELISA. Detects target at ~45 kDa. Cross-reacts with mouse and rat. Ideal for hemolysis, acute phase response, and intestinal permeability studies.

[Learn More](#)

Parameter	Value
Product Name	Haptoglobin (HP) Rabbit mAb
Remarks/Alias	BP; HPA1S; HP2ALPHA2; Haptoglobin (HP)
Species	Human
GeneID (Human)	3240
GeneID	3240
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 307-406 of human Haptoglobin (HP) (P00738)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00738
Protein Weight	45kDa
Shipping	Ice bag