


Anti-Haptoglobin (Hp) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P00738

Item Number: CM0014300



Introduction

Rabbit polyclonal antibody targeting human Haptoglobin (HP), validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Captures free hemoglobin, antioxidant, acute phase protein. Ideal for hemolysis, inflammation, and intestinal permeability research.

[Learn More](#)

Parameter	Value
Product Name	Haptoglobin (HP) Rabbit pAb
Remarks/Alias	BP; HPA1S; HP2ALPHA2; Haptoglobin (HP)
Species	Human
GeneID (Human)	3240
GeneID	3240
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Haptoglobin (NP_005134.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00738
Protein Weight	45kDa
Shipping	Ice bag