

Anti-Has1 Rabbit Polyclonal Antibody For Wb, Elisa - Q92839

Item Number: CM0013405

Introduction

Rabbit polyclonal antibody targeting human Hyaluronan synthase 1 (HAS1). Validated for WB and ELISA, detects a 65 kDa protein. Cross-reacts with human and mouse. Essential for hyaluronan synthesis and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	HAS1 Rabbit pAb
Remarks/Alias	HAS; HAS1
Species	Human
Gene ID	3036
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 80-390 of human HAS1 (NP_001514.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92839
Protein Weight	65kDa
Shipping	Ice bag