

Anti-Has3 Polyclonal Antibody For Wb And Elisa - O00219

Item Number: CM0026800

Introduction

Rabbit polyclonal antibody targeting human HAS3 for WB and ELISA. Cross-reacts with mouse, rat. Hyaluronan synthase 3 essential for hyaluronan synthesis and extracellular matrix remodeling. Suitable for cell adhesion, migration, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HAS3 Rabbit pAb
Remarks/Alias	HAS3
Species	Human
Gene ID (Human)	3038
Gene ID	3038
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-281 of human HAS3 (NP_619515.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00219
Protein Weight	63kDa
Shipping	Ice bag