

## Anti-Has3 Polyclonal Antibody For Wb, If/Icc, Elisa - 000219

Item Number: CM0013919

### Introduction

Rabbit polyclonal antibody against human HAS3, validated for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Suitable for hyaluronan synthesis research.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | HAS3 Rabbit pAb                                                                                                                       |
| Remarks/Alias    | HAS3                                                                                                                                  |
| Species          | Human                                                                                                                                 |
| GeneID (Human)   | 3038                                                                                                                                  |
| GeneID           | 3038                                                                                                                                  |
| Immunogen        | Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 67-281 of human HAS3 (NP_619515.1) |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, IF/ICC, ELISA                                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | O00219                                                                                                                                |
| Protein Weight   | 63kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |