

Anti-Haspin Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tf76

Item Number: CM0020453



Introduction

High-quality rabbit polyclonal antibody against human HASPIN (GSG2). Validated for WB and ELISA, cross-reacts with human and mouse. Research use for cell cycle and chromatin regulation studies.

[Learn More](#)

Parameter	Details
Product Name	HASPIN Rabbit pAb
Remarks/Alias	GSG2; HASPIN
Species	Human
GeneID (Human)	83903
GeneID	83903
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-470 of human HASPIN (NP_114171.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TF76
Protein Weight	88kDa
Shipping	Ice bag