

## Anti-Haus6 Polyclonal Antibody For Wb And Elisa - Q7Z4H7

Item Number: CM0031058



### Introduction

Rabbit polyclonal antibody targeting human HAUS6, a key component of the HAUS augmin complex. Validated for WB and ELISA; cross-reacts with human and mouse. Useful for mitosis and centrosome studies.

[Learn More](#)

| Parameter        | Value                                                                                  |
|------------------|----------------------------------------------------------------------------------------|
| Product Name     | HAUS6 Rabbit pAb                                                                       |
| Remarks/Alias    | Dgt6, FAM29A, HAUS6                                                                    |
| Species          | Human                                                                                  |
| Gene ID          | 54801                                                                                  |
| Immunogen        | Recombinant fusion protein containing amino acids 676-955 of human HAUS6 (NP_060115.3) |
| Source           | Rabbit                                                                                 |
| Category         | Polyclonal Antibodies                                                                  |
| Application      | WB, ELISA                                                                              |
| Cross Reactivity | Human, Mouse                                                                           |
| SWISS            | Q7Z4H7                                                                                 |
| Protein Weight   | 109 kDa                                                                                |
| Shipping         | Ice bag                                                                                |