

Anti-Hax1 Polyclonal Antibody For Wb, If/Icc, Elisa - 000165

Item Number: CM0024855

Introduction

Rabbit polyclonal HAX1 antibody validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Detects HCLS1-associated protein X-1 (32kDa) involved in actin remodeling and cell survival.

[Learn More](#)

Parameter	Value
Product Name	HAX1 Rabbit pAb
Remarks/Alias	SCN3, HS1BP1, HCLSBP1, HAX1
Species	Human
GeneID (Human)	10456
GeneID	10456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-279 of human HAX1 (NP_006109.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00165
Protein Weight	32kDa
Shipping	Ice bag