


Anti-Hb-Egf Rabbit Monoclonal Antibody For Wb, Elisa - Q99075

Item Number: CM0002392



Introduction

Rabbit monoclonal antibody targeting human HB-EGF (DTR, HEGFL) for Western blot and ELISA applications. Validated in human, mouse, and rat samples. Ideal for growth factor and diphtheria toxin receptor studies.

[Learn More](#)

Parameter	Value
Product Name	HB-EGF Rabbit mAb
Remarks/Alias	DTR; DTS; DTSF; HEGFL; HB-EGF
Species	Human
GeneID (Human)	1839
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human HB-EGF (Q99075)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99075
Protein Weight	23kDa
Shipping	Ice bag