


Anti-Hb-Egf Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99075

Item Number: CM0024627



Introduction

Anti-HB-EGF Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Targets proheparin-binding EGF-like growth factor, a key EGFR ligand involved in cardiac development and smooth muscle proliferation. Ideal for growth factor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	HB-EGF Rabbit pAb
Remarks/Alias	DTR, DTS, DTSF, HEGFL, HB-EGF
Species	Human
GeneID (human)	1839
GeneID	1839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-160 of human HB-EGF (NP_001936.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99075
Protein Weight	23kDa
Shipping	Ice bag