


Anti-Hba2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P69905

Item Number: CM0029207



Introduction

This anti-HBA2 rabbit polyclonal antibody is validated for WB, IF-P, and ELISA in human, mouse, and rat. Targets hemoglobin subunit alpha for oxygen transport research.

[Learn More](#)

Parameter	Value
Product Name	HBA2 Rabbit pAb
Remarks/Alias	HBH; ECT7; HBA-T2; HBA2
Species	Human
GeneID (Human)	3040
GeneID	3040
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P69905
Protein Weight	15kDa
Shipping	Ice bag