


Anti-Hbb Polyclonal Antibody For Wb, Elisa - P68871

Item Number: CM0014282



Introduction

Anti-HBB Rabbit Polyclonal Antibody validated for WB and ELISA. Targets hemoglobin subunit beta (P68871) in human and rat. Ideal for oxygen transport and anemia research.

[Learn More](#)

Parameter	Value
Product Name	HBB Rabbit pAb
Remarks/Alias	ECYT6; CD113t-C; beta-globin; HBB
Species	Human
GeneID (human)	3043
GeneID	3043
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 47-147 of human HBB (NP_000509.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P68871
Protein Weight	16kDa
Shipping	Ice bag