


Anti-Hbd (Hemoglobin Subunit Delta) Rabbit Polyclonal Antibody For Wb And Elisa - P02042

Item Number: CM0027890



Introduction

Rabbit polyclonal antibody against human hemoglobin subunit delta (HBD) for Western blot (WB) and ELISA applications. Detects mouse and rat orthologs. Suitable for oxygen transport and hemoglobin research.

[Learn More](#)

Parameter	Value
Product Name	HBD Rabbit pAb
Remarks / Alias	HBD
Species	Human
GeneID (Human)	3045
GeneID	3045
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human HBD (NP_000510.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P02042
Protein Weight	16kDa
Shipping	Ice bag