


Anti-Hb γ 2 Rabbit Polyclonal Antibody For Wb And Elisa - P69892

Item Number: CM0029149



Introduction

High-quality HBG2 rabbit polyclonal antibody validated for Western blot and ELISA, targeting human hemoglobin subunit gamma-2 (fetal hemoglobin). Ideal for research on hemoglobinopathies, oxygen transport, and red blood cell biology.

[Learn More](#)

Parameter	Value
Product Name	HBG2 Rabbit pAb
Remarks/Alias	TNCY; HBG-T1; HBG2
Species	Human
Gene ID (Human)	3048
Gene ID	3048
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HBG2 (NP_000175.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P69892
Protein Weight	16kDa
Shipping	Ice bag