


Anti-Hb γ 1 Rabbit Polyclonal Antibody For Wb, Elisa - P69891

Item Number: CM0017613



Introduction

Rabbit polyclonal antibody targeting human HBG1 (fetal hemoglobin gamma-1 chain). Suitable for WB and ELISA applications. Cross-reacts with mouse and rat. Useful for studying fetal hemoglobin development and oxygen transport.

[Learn More](#)

Parameter	Value
Product Name	HBG1 Rabbit pAb
Remarks/Alias	HBGA; HBGR; HBG-T2; HSGGL1; PRO2979; HBG1
Species	Human
Gene ID (Human)	3047
Gene ID	3047
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HBG1 (NP_000550.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P69891
Protein Weight	16kDa
Shipping	Ice bag