

Anti-Hcar1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bxc0

Item Number: CM0019737



Introduction

Rabbit polyclonal antibody targeting HCAR1 (GPR81), reactive with human, mouse, and rat. Suitable for Western blotting and ELISA applications, with a predicted protein weight of 39 kDa. Ideal for studying lactate signaling and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	HCAR1 Rabbit pAb
Remarks/Alias	HCA1; GPR81; LACR1; FKSG80; GPR104; TAGPCR; TA-GPCR; HCAR1
Species	Human
GeneID (Human)	27198
GeneID	27198
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXC0
Protein Weight	39kDa
Shipping	Ice bag