


Anti-Hcfc1 Rabbit Polyclonal Antibody For Wb, Elisa - P51610

Item Number: CM0020920



Introduction

Rabbit polyclonal antibody targeting HCFC1 (Host cell factor 1), a transcriptional coregulator involved in cell cycle control. Validated for WB and ELISA applications; cross-reacts with human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	HCFC1 Rabbit pAb
Remarks/Alias	CFF; HCF; HCF1; HFC1; MRX3; VCAF; HCF-1; MAHCX; XLID3; PPP1R89; HCFC1
Species	Human
GeneID (Human)	3054
GeneID	3054
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human HCFC1 (NP_005325.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51610
Protein Weight	209kDa
Shipping	Ice bag