

Anti-Hcg Beta Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P0Dn86

Item Number: CM0010649



Introduction

Rabbit monoclonal antibody specific for human hCG beta (CGB), suitable for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Ideal for reproductive biology and pregnancy research.

[Learn More](#)

Parameter	Value
Product Name	hCG beta Rabbit mAb
Remarks/Alias	CGB; HCG; hCGB
Species	Human
Gene ID (Human)	93659
Gene ID	93659
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-165 of human CGB5 (NP_000728.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P0DN86
Protein Weight	18kDa
Shipping	Ice bag