


Anti-Hck Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08631

Item Number: CM0022186



Introduction

High-quality rabbit polyclonal antibody against human HCK (p59Hck/p61Hck). Validated for WB, IF/ICC, and ELISA. Ideal for innate immunity, phagocytosis, and cell signaling research. Sourced from CamelBio.

[Learn More](#)

Parameter	Specification
Product Name	HCK Rabbit pAb
Remarks/Alias	JTK9, AIPCV, p59Hck, p61Hck, HCK
Species	Human
GeneID (human)	3055
GeneID	3055
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human HCK (NP_001165600.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P08631
Protein Weight	60kDa
Shipping	Ice bag