

Anti-Hcls1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P14317

Item Number: CM0024423

Introduction

High-quality rabbit polyclonal antibody against human HCLS1 (HS1, p75), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Ideal for antigen receptor signaling and hematopoietic lineage research.

[Learn More](#)

Parameter	Value
Product Name	HCLS1 Rabbit pAb
Remarks/Alias	HS1; p75; CTTNL; IckBP1; HCLS1
Species	Human
Gene ID (Human)	3059
Gene ID	3059
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human HCLS1 (NP_005326.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P14317
Protein Weight	54kDa
Shipping	Ice bag