


Anti-Hcn1 Rabbit Polyclonal Antibody For Wb, Elisa - O60741

Item Number: CM0022350



Introduction

Rabbit polyclonal antibody targeting HCN1 (BCNG-1) for Western blot and ELISA. Validated in human, mouse, and rat. Explore its role in neuronal hyperpolarization-activated currents and pacemaker activity.

[Learn More](#)

Parameter	Value
Product Name	HCN1 Rabbit pAb
Remarks/Alias	BCNG1, DEE24, HAC-2, BCNG-1, EIEE24, GEFSP10, HCN1
Species	Human
GeneID (Human)	348980
GeneID	348980
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 791-890 of human HCN1 (NP_066550.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60741
Protein Weight	99kDa
Shipping	Ice bag