

Anti-Hcov-229E Spike S1 Monoclonal Antibody For Wb, Elisa - P15423

Item Number: CM0007683



Introduction

High-quality rabbit monoclonal antibody against HCoV-229E Spike S1 glycoprotein, specifically recognizing the receptor-binding domain. Validated for Western blot and ELISA, suitable for virology research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	HCoV-229E Spike S1 Rabbit mAb
Gene ID	918758
Immunogen	Recombinant protein of human HCoV-229E Spike S1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-229E
SWISS	P15423
Protein Weight	129kDa
Shipping	Ice bag