

Anti-Hcov-229E Spike S1 Polyclonal Antibody For Wb, Elisa - P15423

Item Number: CM0020393



Introduction

Rabbit polyclonal antibody against HCoV-229E Spike S1, validated for WB and ELISA. Recognizes a synthetic peptide within amino acids 1074-1173. Suitable for investigating coronavirus attachment and entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	HCoV-229E Spike S1 Rabbit pAb
GeneID	918758
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1074-1173 of coronavirus Spike S1 (NP_073551.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-229E
SWISS	P15423
Protein Weight	129kDa
Shipping	Ice bag