

Anti-Hcov-229E Spike S2 Polyclonal Antibody For Wb, Elisa - P15423

Item Number: CM0018894



Introduction

Rabbit polyclonal antibody targeting the S2 subunit of HCoV-229E spike glycoprotein. Validated for Western blot and ELISA, with specific cross-reactivity against HCoV-229E. Ideal for coronavirus research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	HCoV-229E Spike S2 Rabbit pAb
GeneID	918758
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1074-1173 of coronavirus Spike S2 (NP_073551.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-229E
SWISS	P15423
Protein Weight	129kDa
Shipping	Ice bag