

Anti-Hcov-229E Spike S2 Rabbit Monoclonal Antibody For Wb And Elisa - P15423

Item Number: CM0007334



Introduction

Rabbit monoclonal antibody targeting the S2 subunit of HCoV-229E spike glycoprotein, validated for Western blot and ELISA. Suitable for coronavirus entry and fusion research.

[Learn More](#)

Parameter	Value
Product Name	HCoV-229E Spike S2 Rabbit mAb
Gene ID	918758
Immunogen	Recombinant protein of human HCoV-229E Spike S2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-229E
SWISS	P15423
Protein Weight	129kDa
Shipping	Ice bag