

Anti-Hcov-Hku1 Nucleoprotein Rabbit Polyclonal Antibody For Wb, Elisa - Q5Mqc6

Item Number: CM0023820



Introduction

Rabbit polyclonal antibody targeting HCoV-HKU1 nucleoprotein (N protein), suitable for WB and ELISA. Ideal for studying coronavirus replication, assembly, and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	HCoV-HKU1 Nucleoprotein Rabbit pAb
Species	HCoV-HKU1
GeneID	3200423
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of coronavirus Nucleoprotein (YP_173242.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-HKU1
SWISS	Q5MQC6
Protein Weight	49kDa
Shipping	Ice bag