


Anti-Hcov-Hku1 Spike S1 Polyclonal Antibody For Wb, Elisa - Q5Mqd0

Item Number: CM0020423



Introduction

Rabbit polyclonal antibody specific for HCoV-HKU1 Spike S1 protein (immunogen: peptide within aa 100-200). Validated in WB and ELISA. Recognizes HCoV-HKU1 spike glycoprotein. Ideal for studies of viral entry and host interaction. SWISS Q5MQD0.

[Learn More](#)

Parameter	Value
Product Name	HCoV-HKU1 Spike S1 Rabbit pAb
Species	Human coronavirus HKU1
Gene ID	3200426
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of coronavirus Spike S1 (YP_173238.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-HKU1
SWISS	Q5MQD0
Protein Weight	152 kDa
Shipping	Ice bag