


Anti-Hcov-Hku1 Spike S2 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Mqd0

Item Number: CM0021582



Introduction

High-quality rabbit polyclonal antibody targeting the Spike S2 subunit of HCoV-HKU1. Validated for WB and ELISA, this antibody detects viral fusion protein for diagnostic assay development and research on coronavirus entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	HCoV-HKU1 Spike S2 Rabbit pAb
Remarks/Alias	E2, Peplomer protein
GeneID	3200426
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 800-900 of coronavirus Spike S2 (YP_173238.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-HKU1
SWISS	Q5MQD0
Protein Weight	152kDa
Shipping	Ice bag