


Anti-Hcov-NL63 Spike S2 Polyclonal Antibody For Wb, Elisa - Q6Q1S2

Item Number: CM0020477



Introduction

Rabbit polyclonal antibody targeting HCoV-NL63 Spike S2 protein (Spike glycoprotein, E2). Suitable for WB and ELISA applications. Recognizes HCoV-NL63 spike. For coronavirus research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	HCoV-NL63 Spike S2 Rabbit pAb
Gene ID	2943499
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1200-1300 of coronavirus Spike S2 (YP_003767.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-NL63
SWISS	Q6Q1S2
Protein Weight	150kDa
Shipping	Ice bag