

Anti-Hcov-OC43 Nucleoprotein Polyclonal Antibody For Wb, Elisa - P33469

Item Number: CM0020571



Introduction

Rabbit polyclonal antibody targeting HCoV-OC43 Nucleoprotein (N protein) for WB and ELISA. Cross-reacts with HCoV-OC43; immunogen from amino acids 100-200. Suitable for coronavirus nucleocapsid and assembly studies.

[Learn More](#)

Parameter	Value
Product Name	HCoV-OC43 Nucleoprotein Rabbit pAb
Gene ID	39105221
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of coronavirus Nucleoprotein (YP_009555245.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-OC43
SWISS	P33469
Protein Weight	49kDa
Shipping	Ice bag