


Anti-Hcov-OC43 Spike S2 Rabbit Polyclonal Antibody For Wb, Elisa - P36334

Item Number: CM0019775



Introduction

Rabbit polyclonal antibody targeting the S2 subunit of HCoV-OC43 spike glycoprotein (Swiss-Prot P36334). Suitable for Western blot and ELISA applications, validated to detect HCoV-OC43 cross-reactivity. Ideal for coronavirus research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	HCoV-OC43 Spike S2 Rabbit pAb
GeneID	39105218
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 782-881 of coronavirus Spike S2 (YP_009555241.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	HCoV-OC43
SWISS	P36334
Protein Weight	150kDa
Shipping	Ice bag