

Anti-Hcrt Rabbit Polyclonal Antibody For Wb, Elisa - 043612

Item Number: CM0027092



Introduction

Rabbit polyclonal antibody against human HCRT (Orexin precursor), suitable for Western blot and ELISA. Validated in human, mouse, and rat species. Ideal for studying sleep-wake regulation and energy metabolism.

[Learn More](#)

Parameter	Value
Product Name	HCRT Rabbit pAb
Remarks/Alias	OX; PPOX; NRCLP1; HCRT
Species	Human
Gene ID (Human)	3060
Gene ID	3060
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 34-66 of human HCRT (NP_001515.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43612
Protein Weight	13kDa
Shipping	Ice bag