


Anti-Hcrtr2 Polyclonal Antibody For Wb, Elisa - 043614

Item Number: CM0028266



Introduction

High-quality rabbit polyclonal antibody targeting the human orexin receptor type 2 (HCRT2). Validated for Western blot and ELISA, cross-reacts with human and mouse. Suitable for studying sleep-wake regulation and glucose homeostasis.

[Learn More](#)

Parameter	Specification
Product Name	HCRT2 Rabbit pAb
Remarks/Alias	OX2R; OXR2; ORXR2; HCRT2
Species	Human
GeneID (Human)	3062
GeneID	3062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-444 of human HCRT2 (NP_001517.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43614
Protein Weight	51kDa
Shipping	Ice bag