

Anti-Hd1 Polyclonal Antibody For Wb And Elisa - Q9Fdx8

Item Number: CM0020410



Introduction

Anti-Hd1 rabbit polyclonal antibody validated for Western blot and ELISA applications in *Oryza sativa*. Detects Zinc finger protein HD1, a transcription factor regulating flowering. Suitable for plant developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Hd1 Rabbit pAb
Remarks/Alias	Hd1; SE1; OsHd1; Hd1
Species	<i>Oryza sativa</i>
Gene ID	4340746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of <i>Oryza sativa</i> Hd1 (Q9FDX8)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q9FDX8
Protein Weight	43kDa
Shipping	Ice bag