

Anti-Hdac1 Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q13547

Item Number: CM0013561



Introduction

KD-validated rabbit polyclonal antibody targeting HDAC1. Suitable for WB, IP, IF/ICC, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. A reliable tool for epigenetics and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] HDAC1 Rabbit pAb
Remarks/Alias	HD1; RPD3; KDAC1; GON-10; RPD3L1; C1
Species	Human
GeneID (Human)	3065
GeneID	3065
Immunogen	Synthetic peptide Recombinant fusion protein containing a sequence corresponding to amino acids 393-482 of human HDAC1 (NP_004955.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13547
Protein Weight	55 kDa
Shipping	Ice bag