

Anti-Hdac1 Rabbit Monoclonal Antibody For Ihc-P, If/Icc, Elisa - Q13547

Item Number: CM0002414



Introduction

Anti-HDAC1 Rabbit mAb (55kDa) validated for IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat HDAC1. Ideal for studying chromatin remodeling, transcription regulation, and epigenetic pathways.

[Learn More](#)

Parameter	Value
Product Name	HDAC1 Rabbit mAb
Remarks/Alias	HD1; RPD3; KDAC1; GON-10; RPD3L1; C1
Species	Human
Gene ID	3065
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human HDAC1 (Q13547)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13547
Protein Weight	55kDa
Shipping	Ice bag