

Anti-Hdac11 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Db2

Item Number: CM0024959

Introduction

High-quality rabbit polyclonal antibody against human HDAC11 (SWISS: Q96DB2) validated for WB and ELISA. Ideal for epigenetic research on histone deacetylation and transcriptional repression.

[Learn More](#)

Parameter	Value
Product Name	HDAC11 Rabbit pAb
Remarks/Alias	HD11; HDAC11
Species	Human
GeneID (Human)	79885
GeneID	79885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-347 of human HDAC11 (NP_079103.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96DB2
Protein Weight	39kDa
Shipping	Ice bag