


Anti-Hdac2 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - Q92769

Item Number: CM0006601



Introduction

Rabbit monoclonal antibody targeting human HDAC2 (histone deacetylase 2). Suitable for WB, IHC-P, IF/ICC, IP, ELISA, and ChIP. Cross-reacts with mouse and rat. Ideal for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	HDAC2 Rabbit mAb
Remarks/Alias	HD2; RPD3; YAF1; KDAC2; HDAC2
Species	Human
Gene ID	3066
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 389-488 of human HDAC2 (Q92769).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92769
Protein Weight	55kDa
Shipping	Ice bag