


Anti-Hdac2 Rabbit Monoclonal Antibody For Wb, Elisa - Q92769

Item Number: CM0003936



Introduction

Recombinant rabbit monoclonal antibody against human HDAC2, suitable for Western blot and ELISA. Cross-reacts with human, mouse, rat. Ideal for epigenetic research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	HDAC2 Rabbit mAb
Remarks/Alias	HD2, RPD3, YAF1, KDAC2, HDAC2
Species	Human
Gene ID (Human)	3066
Gene ID	3066
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HDAC2 (NP_001518.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92769
Protein Weight	55kDa
Shipping	Ice bag