


Anti-Hdac3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15379

Item Number: CM0013913



Introduction

Rabbit polyclonal antibody targeting HDAC3, suited for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat HDAC3. Ideal for chromatin biology and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	HDAC3 Rabbit pAb
Remarks/Alias	HD3; RPD3; KDAC3; RPD3-2; HDAC3
Species	Human
Gene ID (Human)	8841
Gene ID	8841
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 329-428 of human HDAC3 (NP_003874.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15379
Protein Weight	49 kDa
Shipping	Ice bag