


## Anti-Hdac3 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - O15379

Item Number: CM0009519



### Introduction

High-quality Anti-HDAC3 Rabbit Monoclonal Antibody validated for WB, IF/ICC, IHC-P, ELISA with reactivity across human, mouse, rat, and monkey. Ideal for epigenetic, chromatin biology, and transcriptional regulation studies.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | HDAC3 Rabbit mAb                                                                                                   |
| Remarks/Alias    | HD3; RPD3; KDAC3; RPD3-2                                                                                           |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 8841                                                                                                               |
| GeneID           | 8841                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 299-428 of human HDAC3 (NP_003874.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                              |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat, Monkey                                                                                          |
| SWISS            | O15379                                                                                                             |
| Protein Weight   | 49 kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |