


Anti-Hdac3 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - O15379

Item Number: CM0003576



Introduction

Rabbit monoclonal HDAC3 antibody validated for Western blot, immunoprecipitation, and ELISA with cross-reactivity to human, mouse, and rat. The 49 kDa histone deacetylase 3 (HDAC3) is essential for chromatin remodeling and gene silencing. Suitable for epigenetic and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] HDAC3 Rabbit mAb
Remarks/Alias	HD3; RPD3; KDAC3; RPD3-2; HDAC3
Species	Human
GeneID (Human)	8841
GeneID	8841
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 329-428 of human HDAC3 (NP_003874.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15379
Protein Weight	49kDa
Shipping	Ice bag