

Anti-Hdac4 Polyclonal Antibody For Wb, Elisa - P56524

Item Number: CM0012815



Introduction

Rabbit polyclonal antibody against human HDAC4, raised with a synthetic peptide. Validated for WB and ELISA, this antibody detects the 119kDa histone deacetylase 4, a key epigenetic regulator.

[Learn More](#)

Parameter	Value
Product Name	HDAC4 Rabbit pAb
Remarks/Alias	HD4; AHO3; BDMR; HDACA; HA6116; HDAC-4; HDAC-A; NEDCHF; NEDCHID; HDAC4
Species	Human
Gene ID (Human)	9759
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 530-630 of human HDAC4 (NP_006028.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P56524
Protein Weight	119kDa
Shipping	Ice bag