

Anti-Hdac4 Polyclonal Antibody For Wb, Ihc-P, Elisa - P56524

Item Number: CM0013964



Introduction

Rabbit polyclonal antibody against human HDAC4 (SWISS: P56524). Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for epigenetic and cancer research. Protein weight 119 kDa.

[Learn More](#)

Parameter	Value
Product Name	HDAC4 Rabbit pAb
Remarks/Alias	HD4; AHO3; BDMR; HDACA; HA6116; HDAC-4; HDAC-A; NEDCHF; NEDCHID; HDAC4
Species	Human
Gene ID (Human)	9759
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 451-650 of human HDAC4 (NP_006028.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56524
Protein Weight	119kDa
Shipping	Ice bag