

Anti-Hdac5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uql6

Item Number: CM0015733



Introduction

This HDAC5 rabbit polyclonal antibody is validated for WB, IHC-P, IF/ICC, and ELISA applications. It detects HDAC5 in human, mouse, and rat samples, making it suitable for chromatin regulation and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	HDAC5 Rabbit pAb
Remarks/Alias	HD5; NY-CO-9; HDAC5
Species	Human
Gene ID (Human)	10014
Gene ID	10014
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 500-600 of human HDAC5 (NP_005465.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQL6
Protein Weight	122kDa
Shipping	Ice bag